

Product Datasheet

NEFH antibody (orb388517)

Description

Mouse monoclonal antibody to NEFH

Species/Host

Mouse

Reactivity

Bovine, Feline, Gallus, Guinea pig, Human, Mouse, Porcine, Rabbit, Rat, Rodent

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunohistochemical
staining of human
Ce...

Immunogen

Crude neurofilament preparation from porcine spinal cord

Target

NEFH

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Neurofilament, phospho (NF-H) (Neuronal Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

200kDa

Uniprot ID

P12036

Entrez



4744

Dilution Range

FACS: 0.5-1 µg/million cells, IHC-P: 0.25-0.5 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical
staining of rat Cere...