


Product Datasheet

NEFL antibody (orb388520)

Description	Mouse monoclonal antibody to NEFL
Species/Host	Mouse
Reactivity	Bovine, Gallus, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	Recombinant human NEFL protein
Target	NEFL
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Neurofilament (NF-L) (Neuronal Marker)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	68kDa
Uniprot ID	P07196
Entrez	4747
Dilution Range	IF/ICC: 1-2 µg/ml, FACS: 0.5-1 µg/million cells, IHC-P: 0.25-0.5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat Cere...