

Product Datasheet

NGFR antibody (orb388521)

Description

Mouse monoclonal antibody to NGFR


Species/Host

Mouse

Reactivity

Feline, Ferret, Human, Monkey, Primate, Rabbit

Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications

Immunohistochemical
staining of human
Me...

Immunogen

NGFR from A875 melanoma cells

Target

NGFR

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

NGF-Receptor (p75) / CD271 (Soft Tissue Tumor Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

75kDa

Uniprot ID

P08138

Entrez

4804

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.