

Product Datasheet

PIP antibody (orb388557)

Description

Mouse monoclonal antibody to PIP

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications

Immunogen

Recombinant human GCDFP-15 protein fragment (around aa 41-146) (exact sequence is proprietary)

Target

PIP

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

GCDFP-15 (Gross Cystic Disease Fluid Protein 15) (Breast Marker)

Isotype

Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

15kDa

Uniprot ID

[P12273](#)

Entrez


[5304](#)

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.



Western blot analysis of human Pancreas ...



Immunohistochemical staining of human Br...