

## Product Datasheet

### CYCS antibody (orb388562)

**Description**

Mouse monoclonal antibody to CYCS

**Species/Host**

Mouse

**Reactivity**

Aves, Canine, Drosophila, Equine, Frog, Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested**

FC, IHC-P, WB

**Applications**

Immunohistochemical  
staining of human  
He...

**Immunogen**

Recombinant cytochrome c protein

**Target**

CYCS

**Preservatives**

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

**Concentration**

Purified Ab with BSA and Azide at 200ug/ml

**Storage**

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

**Note**

For research use only

**Application notes**

Cytochrome C (Mitochondrial Marker)

**Isotype**

Mouse / IgG2b, kappa

**Clonality**

Monoclonal

**MW**

15kDa

**Uniprot ID**

**P99999**

**Entrez**

**54205**

**Dilution Range**

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 0.25-0.5 µg/ml

**Expiration Date**

12 months from date of receipt.

