

Product Datasheet

TMEM16A antibody (orb388576)

Description

Mouse monoclonal antibody to TMEM16A

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical
staining of human
Gl...

Immunogen

A synthetic peptide from human DOG-1 protein (aa 904-986), conjugated to a carrier protein. (

Target

TMEM16A

Preservatives

200ug/ml of Ab purified from Bioreactor
Concentrate by Protein A/G. Prepared in 10mM PBS
with 0.05% rAlbumin & 0.05% azide. Also available
WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8 °C. Antibody
without azide - store at -20 to -80 °C. Antibody is
stable for 24 months. Non-hazardous. No MSDS
required.

Note

For research use only

Application notes

DOG-1 / TMEM16A / ANO1 (Gastrointestinal Stromal
Tumor Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

~114kDa

Uniprot ID

[Q5XXA6](#)

Entrez

[55107](#)

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml,
IHC-P: 0.25-0.5 µg/ml

Expiration Date

12 months from date of receipt.