

Product Datasheet

ACPP antibody (orb388584)

Description

Mouse monoclonal antibody to ACPP

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant full-length human ACPP protein
Target	ACPP
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Prostate Specific Acid Phosphatase (PSAP)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	52kDa
Uniprot ID	P15309
Entrez	55
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-Fr: 0.5-1.0 µg/ml
Expiration Date	12 months from date of receipt.