

Product Datasheet

BAX antibody (orb388598)

Description

Mouse monoclonal antibody to BAX

Species/Host

Mouse

Reactivity


Human, Monkey

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications


Immunohistochemical staining of Hodgkin'...

Immunogen

A synthetic peptide, aa 3-16 (Cys-GSGEQPRGGGPTSS) of human bax protein.

Target

BAX

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Bax (Apoptosis Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

21kDa

Uniprot ID

[Q07812](#)

Entrez

[581](#)

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.