

Product Datasheet

BCL2L1 antibody (orb388623)

Description

Mouse monoclonal antibody to BCL2L1

Species/Host

Mouse

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

A synthetic peptide, aa 3-14 (Cys-QSNRELVVDFLS) of human Bcl-X protein (BX006 & 2H12)

Target

BCL2L1

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Bcl-X (Apoptosis Marker)

Isotype

Mouse / IgG2a, kappa + Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

27kDa

Uniprot ID

[Q07817](#)

Entrez


[598](#)

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.



Western blot analysis of Jurkat Cell Lys...



Immunohistochemical staining of human Ho...