


Product Datasheet

BCL6 antibody (orb388624)

Description	Mouse monoclonal antibody to BCL6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human bcl-6 protein fragment (around aa256-389) (exact sequence is proprietary)
Target	BCL6
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Bcl-6 (Follicular Lymphoma Marker)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	95kDa
Uniprot ID	P41182
Entrez	604
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, WB: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 Cell Lysa...