

Product Datasheet

BCL6 antibody (orb388625)

Description

Mouse monoclonal antibody to BCL6

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical
staining of human
To...

Immunogen

Recombinant human bcl-6 protein fragment (aa256-389) (Exact sequence is proprietary)

Target

BCL6

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Bcl-6 (Follicular Lymphoma Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

95kDa

Uniprot ID

[P41182](#)

Entrez

[604](#)

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, IHC-P: 1-2 µg/ml

Expiration Date

12 months from date of receipt.