


Product Datasheet

S100A1 antibody (orb388627)

Description	Mouse monoclonal antibody to S100A1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	Purified human full-length native S-100 protein conjugated to methylated BSA
Target	S100A1
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	S100A (Melanoma Marker)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	10-12kDa
Uniprot ID	P23297
Entrez	6271
Dilution Range	IF/ICC: 1-2 µg/ml, IHC-P1-2 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Me...