


Product Datasheet

S100B antibody (orb388641)

Description	Rabbit monoclonal antibody to S100B
Species/Host	Rabbit
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full-length human S100B protein
Target	S100B
Preservatives	200ug/ml of Ab purified by Protein A. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide store at 2 to 8 °C. Antibody without azide store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	S100B (Astrocyte and Melanoma Marker)
Isotype	Rabbit / IgG, kappa
Clonality	Monoclonal
MW	10-12kDa
Uniprot ID	P04271
Entrez	6285
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 0.25-0.5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Me...