

Product Datasheet

FBXW2 antibody (Biotin) (orb38865)

Description	FBXW2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human F-box/WD repeat-containing protein 2 protein (72-340AA)
Target	FBXW2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UKT8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.