


Product Datasheet

BSG antibody (orb388684)

Description	Mouse monoclonal antibody to BSG
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Immunogen	Human T cell line Molt 13
Target	BSG
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	CD147 / EMMPRIN / Neurothelin
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	35kDa (non-reduced); 40kDa (reduced)
Uniprot ID	P35613
Entrez	682
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 1:100-1:200, IHC-P: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Re...