

Product Datasheet

TGFA antibody (orb388702)

Description	Mouse monoclonal antibody to TGFA
Species/Host	Mouse
Reactivity	Human, Rabbit, Zebrafish
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Immunogen	A 10-amino acid synthetic peptide (aa34-43; PPVAAAVVSH) from human TGF α .
Target	TGFA
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	TGF-alpha (Transforming Growth Factor alpha)
Isotype	Mouse / IgM, kappa
Clonality	Monoclonal
MW	~6kDa
Uniprot ID	P01135
Entrez	7039
Dilution Range	FACS: 0.5-1 μ g/million cells, IF/ICC: 1-2 μ g/ml, IHC-P: 2-4 μ g/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Te...