



Product Datasheet

TGM2 antibody (orb388709)


Description	Mouse monoclonal antibody to TGM2
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full-length human TGM2 protein
Target	TGM2
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Transglutaminase II (TGM2)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	77-85kDa
Uniprot ID	P21980
Entrez	7052
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human To...



Immunohistochemical staining of human Pl...



Immunohistochemical staining of human Ut...