

Product Datasheet

TIMP2 antibody (orb388716)

Description

Mouse monoclonal antibody to TIMP2

Species/Host

Mouse

Reactivity

Bovine, Guinea pig, Human, Mouse, Rabbit, Rat, Zebrafish

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications
Immunogen

A synthetic peptide from the N-terminus of human TIMP2 protein (aa56-66) (DSGNDIYGNPI)

Target

TIMP2

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

TIMP2 (Tissue Inhibitor of Metalloproteinase 2)

Isotype

Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

21kDa

Uniprot ID

P16035

Entrez


7077

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 1-2 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of human Pl...