

Product Datasheet

TNF antibody (orb388733)

Description	Mouse monoclonal antibody to TNF
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, NeA
Immunogen	Recombinant human TNF-alpha
Target	TNF
Preservatives	200ug/ml of Ab purified by Protein L. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	TNF-alpha (Tumor Necrosis Factor alpha)
Isotype	Mouse / IgM, kappa
Clonality	Monoclonal
MW	17kDa
Uniprot ID	P01375
Entrez	7124
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml
Expiration Date	12 months from date of receipt.