

Product Datasheet

TYRP1 antibody (orb388763)

Description

Mouse monoclonal antibody to TYRP1

Species/Host

Mouse

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical
staining of human
Me...

Immunogen

Recombinant full-length human TYRP1 protein

Target

TYRP1

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Tyrosinase-Related Protein-1 (TYRP-1) (Melanoma Marker)

Isotype

Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

75kDa

Uniprot ID

P17643

Entrez

7306

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, IHC-P: 3-6 µg/ml

Expiration Date

12 months from date of receipt.