



Product Datasheet

VIM antibody (orb388775)


Description	Mouse monoclonal antibody to VIM
Species/Host	Mouse
Reactivity	Bovine, Canine, Feline, Gallus, Goat, Human, Porcine
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant full-length human vimentin protein
Target	VIM
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Vimentin (Mesenchymal Cell Marker)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	57-60kDa
Uniprot ID	P08670
Entrez	7431
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, WB: 0.25-0.5 µg/ml, IHC-P: 0.25-0.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Raji Cell Lysat...



Immunohistochemical staining of human Me...



Flow cytometric analysis of Jurkat Cells...