

Product Datasheet

WT1 antibody (orb388783)

Description

Mouse monoclonal antibody to WT1

Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

Recombinant full length human WT1 protein (WT1/857); Recombinant fragment aa1-181 of human WT1 (6F-H2)

Target

WT1

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Wilm's Tumor 1 (WT1) (Wilm's Tumor & Mesothelial Marker)

Isotype

Mouse / IgG1, kappa + Mouse / IgG1, kappa

Clonality

Monoclonal

MW

47-55kDa

Uniprot ID

P19544

Entrez


7490

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human Fe...



Immunohistochemical staining of mouse Ki...



Immunohistochemical staining of rat Kidn...