

Product Datasheet

ZAP70 antibody (orb388789)

Description

Mouse monoclonal antibody to ZAP70

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

Recombinant full-length human ZAP70 protein (ZAP70/528); Recombinant ZAP-70 protein including residues 1-254 and encompassing SH2 domains of human ZAP70 (2F3.2)

Target

ZAP70

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

ZAP70 (Chronic Lymphocytic Leukemia Marker)

Isotype

Mouse / IgG2a, kappa + Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

70kDa

Uniprot ID

[P43403](#)

Entrez



[7535](#)

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of human To...