


Product Datasheet

PAX8 antibody (orb388795)

Description	Mouse monoclonal antibody to PAX8
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant fragment (around aa 60-261) of human PAX8 protein (exact sequence is proprietary)
Target	PAX8
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	PAX8 (Renal Cell Marker)
Isotype	Mouse / IgG2b, kappa
Clonality	Monoclonal
MW	62kDa
Uniprot ID	Q06710
Entrez	7849
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, WB: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Raji Cell Lysat...