

Product Datasheet

MAML2 antibody (orb388805)

Description

Mouse monoclonal antibody to MAML2

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications

Immunogen

Recombinant human MAML2 fragment

Target

MAML2

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

MAML2 (Mastermind Like Transcriptional Coactivator 2)

Isotype

Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

125kDa

Uniprot ID

Q8IZL2

Entrez


84441

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 0.5-1.0 µg/ml

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human Pa...



Immunohistochemical staining of human Co...



Immunohistochemical staining of human Bl...