

Product Datasheet

CD86 antibody (orb388890)

Description

Mouse monoclonal antibody to CD86

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, Functional Studies, IF, IHC-P

Applications
Immunogen

ARH-77 (B-lymphoblastoid cell line)

Target

CD86

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

CD86 (Dendritic Cells Maturation Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

70kDa

Uniprot ID

[P42081](#)

Entrez


[942](#)

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1 µg/ml, IP: 0.5-1 µg/500 µg protein Lysate, IHC-P: 2-4 µg/ml

Expiration Date

12 months from date of receipt.



Flow cytometric analysis of human PBMC's...



Immunohistochemical staining of human Es...