

Product Datasheet

TNFRSF8 antibody (orb388893)

Description

Mouse monoclonal antibody to TNFRSF8

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical
staining of human
Ho...

Immunogen

Human CD30 recombinant protein

Target

TNFRSF8

Preservatives


200ug/ml of Ab Purified from Bioreactor
Concentrate by Protein A/G. Prepared in 10mM PBS
with 0.05% rAlbumin & 0.05% azide. Also available
WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody
without azide - store at -20 to -80°C. Antibody is
stable for 24 months. Non-hazardous. No MSDS
required.



Immunohistochemical
staining of human
Ho...

Note

For research use only

Application notes

CD30 / TNFRSF8 (Hodgkin & Reed-Sternberg Cell
Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

105-120kDa

Uniprot ID

P28908

Entrez

943

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, IHC-
P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.