

Product Datasheet

NAPSA antibody (orb388903)

Description

Mouse monoclonal antibody to NAPSA

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications
Immunogen

Recombinant human Napsin-A protein fragment (around aa189-299) (exact sequence is proprietary)

Target

NAPSA

Preservatives

200ug/ml of Ab Purified from rabbit anti-serum by Protein A. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Napsin A (Lung Adenocarcinoma Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

37kDa

Uniprot ID

[O96009](#)

Entrez


[9476](#)

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, WB: 1-2 µg/ml, IHC-P: 1-2 µg/ml

Expiration Date

12 months from date of receipt.



Western blot analysis of K562, HEK293 an...



Immunohistochemical staining of human Lu...