



Product Datasheet

CD63 antibody (orb388929)


Description	Mouse monoclonal antibody to CD63
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Immunogen	Smooth plasma membrane fraction of MeWo cells
Target	CD63
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	CD63 (Late Endosomes Marker)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	26kDa (core protein); 30-60kDa (glycosylated)
Uniprot ID	P08962
Entrez	967
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human MeWo cells



Immunohistochemical staining of mouse spleen tissue



Western blot analysis of human Spleen Lysate