

Product Datasheet

Phosphotyrosine antibody (orb388991)

Description Mouse monoclonal antibody to Phosphotyrosine

Species/Host Mouse

Reactivity Other

Conjugation Unconjugated

Tested Applications FC, IF, IHC, IP, WB

Immunogen Phosphotyrosine conjugated to KLH

Preservatives 200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration Purified Ab with BSA and Azide at 200ug/ml

Storage Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note For research use only

Application notes Phosphotyrosine (P-Tyr)

Isotype Mouse / IgG2b

Clonality Monoclonal


MW Depends upon the phosphorylated target

Uniprot ID Not Applicable

Entrez Not Applicable

Dilution Range FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, WB: 0.5-1 µg/ml, IP: 1-2 µg/500 µg protein Lysate, IHC-P: 1-2 µg/ml

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human Br...