


Product Datasheet

KRT5 antibody (orb388998)

Description	Mouse monoclonal antibody to KRT5
Species/Host	Mouse
Reactivity	Bovine, Canine, Ferret, Hamster, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Cytoskeletal preparation from HeLa cells
Target	KRT5 & KRT8
Preservatives	0.5 ml at 200 µg/ml of antibody Purified from Bioreactor Concentrate by Protein A/G chromatography. Supplied in 10mM PBS, pH 7.4 with 0.05% rAlbumin and 0.05% sodium azide. Contact us if you require this antibody in a different format.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Store at 2 to 8°C. Antibody is stable for 24 months.
Note	For research use only
Application notes	Cytokeratin 5/8
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	58kDa (CK5); 52.5kDa (CK8)
Uniprot ID	P05787 , P13647
Entrez	3856 , 3852
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1.0 µg/ml, WB: 0.5-1.0 µg/ml, IHC-P: 0.5-1.0 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Co...