


Product Datasheet

CA19-9 antibody (orb389027)

Description	Mouse monoclonal antibody to CA19-9
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Immunogen	Precipitin lines obtained after immuno-diffusion using MAb 116-NS-19-9 and mucins isolated from an ovarian cyst of a Lewis A+B- patient (0Le).
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	CA19-9 / Sialyl Lewisa (GI Tumor Marker)
Isotype	Mouse / IgM, kappa
Clonality	Monoclonal
MW	> 400kDa
Uniprot ID	Not Known
Entrez	Not Known
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1.0 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Co...