


Product Datasheet

Cytokeratin antibody (orb389029)

Description

Mouse monoclonal antibody to Cytokeratin


Species/Host

Mouse

Reactivity

Bovine, Frog, Goat, Guinea pig, Human, Marmoset, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications

Immunohistochemical
staining of human
Co...

Immunogen

Keratin-enriched preparation from cultured human A431

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Cytokeratin, Multi (Epithelial Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

Multiple

Uniprot ID

[P04259](#), [P19013](#), [P02538](#), [P48668](#), [P13647](#)

Entrez

[3860](#), [3875](#), [3851](#), [3852](#), [3856](#), [286887](#), [3858](#), [3853](#), [3854](#)

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1.0 µg/ml, WB: 0.5-1.0 µg/ml, IHC-P: 0.5-1.0 µg/ml

Expiration Date

12 months from date of receipt.