

Product Datasheet

Progesterone antibody (orb389050)

Description

Mouse monoclonal antibody to Progesterone

Species/Host

Mouse

Reactivity

All

Conjugation

Unconjugated

Tested

ELISA, ICC, RIA

Applications
Immunogen

Progesterone- 11a-hemisuccinate conjugated to bovine serum albumin (Pr11a-HMS-BSA)

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C.
Antibody without azide - store at -20 to -80°C.
Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Progesterone

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

314.46Da

Uniprot ID

Not Applicable

Entrez

Not Applicable

Dilution Range

ELISA (For coating, order Ab without BSA);
Immunohistochemistry (Not assessed);Radioimmunoassay (RIA)

Expiration Date

12 months from date of receipt.