

Product Datasheet

Macrophage & Histiocytoma Marker antibody (orb389052)

Description

Mouse monoclonal antibody to Macrophage


Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested


IHC-P

Applications

Immunohistochemical
staining of human
Hi...

Immunogen

Membrane preparation from human
hepatocytes


Preservatives

200ug/ml of Ab Purified from Bioreactor
Concentrate by Protein A/G. Prepared in
10mM PBS with 0.05% azide. Also available
WITHOUT azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C.
Antibody without azide - store at -20 to -
80°C. Antibody is stable for 24 months.
Non-hazardous. No MSDS required.

Immunohistochemical
staining of human
Hi...

Note

For research use only

Application notes

Macrophage & Histiocytoma Marker

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

125kDa -135kDa

Uniprot ID

Not Known

Entrez

Not Known

Dilution Range

IHC-P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.