

Product Datasheet

Multiple antibody (orb389081)

Description

Mouse monoclonal antibody to Multiple

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications
Immunogen

Solubilized keratin extract from human stratum corneum

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Cytokeratin, Basic (Type II or HMW) (Epithelial Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

67kDa (CK1); 58kDa (CK5); 56.5kDa (CK10); 50kDa (CK14)

Uniprot ID

[P04264](#), [P02533](#), [P13645](#), [P13647](#)

Entrez



[3848](#), [3858](#), [3852](#), [3861](#)

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1.0 µg/ml, WB: 0.5-1.0 µg/ml, IHC-P: 0.25-0.5 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of human Pr...