



Product Datasheet

Mitochondria antibody (orb389091)

Description	Mouse monoclonal antibody to Mitochondria
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	Semi-purified mitochondrial preparation
Target	Not Known
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Mitochondria (Marker for Human Cells)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	60kDa
Uniprot ID	Not Known
Entrez	Not Known
Dilution Range	IF/ICC: 0.5-1 µg/ml, WB: 0.25-0.5 µg/ml, IHC-P: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human
To...



Immunohistochemical staining of human
Pa...