

Product Datasheet

KRT77 antibody (orb389098)

Description

Mouse monoclonal antibody to KRT77

Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

Recombinant human KRT77 protein

Target

KRT77

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Cytokeratin, Acidic (Type I or LMW) (Epithelial Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

56.5kDa (CK10); 50kDa (CK14); 50kDa (CK15); 48kDa (CK16); 40kDa (CK19)

Uniprot ID

[Q7Z794](#)

Entrez


[374454](#)

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human Sk...



Immunohistochemical staining of rat Ovid...



Immunohistochemical staining of rat Stom...