

Product Datasheet

N Cadherin antibody (orb389106)

Description

Mouse monoclonal antibody to N Cadherin

Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Recombinant human N-cadherin cytoplasmic domain
Target	CDH2
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	N-Cadherin / Cadherin-2 / CD325 (NCAD)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	130-140kDa
Uniprot ID	P19022
Entrez	1000
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, WB: 0.5-1.0 µg/ml
Expiration Date	12 months from date of receipt.