

Product Datasheet

GPC3 antibody (orb389119)

Description

Mouse monoclonal antibody to GPC3

Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications
Immunogen

Recombinant full-length human GPC3 protein

Target

GPC3

Preservatives

200ug/ml of Ab purified by Protein A Column.
Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Glypican-3 (GPC3) (Hepatocellular Carcinoma Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW


67kDa

Uniprot ID
P51654
Entrez
2719
Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of human He...