


Product Datasheet

GPC3 antibody (orb389120)

Description	Rabbit monoclonal antibody to GPC3
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Immunogen	Recombinant human full-length GPC3 protein
Target	GPC3
Preservatives	200ug/ml of Ab Purified by Protein A. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Glypican-3 (GPC3) (Hepatocellular Carcinoma Marker)
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	67kDa
Uniprot ID	P51654
Entrez	2719
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Fe...