

Product Datasheet

pos5 antibody (FITC) (orb38919)

Description	pos5 antibody (FITC)
Species/Host	Rabbit
Reactivity	Fungi
Conjugation	FITC
Immunogen	Recombinant Schizosaccharomyces pombe NADH kinase pos5, mitochondrial protein (43-386AA)
Target	pos5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UT98
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.