

Product Datasheet

PTH antibody (orb389203)

Description

Mouse monoclonal antibody to PTH

Species/Host

Mouse

Reactivity

Bovine, Canine, Deer, Human, Porcine, Primate, Rabbit, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

A synthetic peptide around amino acids 1-34 of human mature-PTH-polypeptide (exact sequence is proprietary)

Target

PTH

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Parathyroid Hormone (PTH) (N-Terminal)

Isotype

Mouse / IgG2b, kappa

Clonality

Monoclonal

MW

9kDa

Uniprot ID

[P01270](#)

Entrez


[5741](#)

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1.0 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of human Pa...