

Product Datasheet

Calprotectin antibody (orb389243)

Description


Mouse monoclonal antibody to Calprotectin


Species/Host

Mouse

Reactivity

Bovine, Canine, Equine, Feline, Goat, Guinea pig, Human, Monkey, Mouse, Porcine, Primate, Rabbit, Rat


Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications

Immunohistochemical
staining of human
To...

Immunogen

Affinity Purified monocyte membrane preparation

Target

S100A8

Preservatives

200ug/ml of Ab Purified from Bioreactor
Concentrate by Protein A/G. Prepared in 10mM
PBS with 0.05% rAlbumin & 0.05% azide. Also
available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody
without azide - store at -20 to -80°C. Antibody is
stable for 24 months. Non-hazardous. No MSDS
required.

Note

For research use only

Application notes

Macrophage L1 Protein

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

MAC387

Uniprot ID

P06702, P05109

Entrez

6280, 6279

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml,
IHC-P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.