



Product Datasheet

PGP9.5 antibody (orb389260)

Description	Mouse monoclonal antibody to PGP9.5
Species/Host	Mouse
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Native Uchl1 (PGP9.5) protein from brain
Target	Uchl1
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	PGP9.5 / Uchl1 (pan-Neuronal Marker)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	20-30kDa
Uniprot ID	P09936
Entrez	7345
Dilution Range	WB: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Ce...



Western blot analysis of human