

Product Datasheet

vWF antibody (orb389272)

Description

Mouse monoclonal antibody to vWF

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC, IP, WB

Applications
Immunogen

Recombinant human vWF fragment spanning around aa 845-949 (exact sequence is proprietary)

Target

VWF

Preservatives

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

von Willebrand Factor / Factor VIII Related-Ag (Endothelial Marker)

Isotype

Mouse / IgG1, kappa (3E2D10) + Mouse / IgG1, kappa (VWF635)

Clonality

Monoclonal

MW

250kDa

Uniprot ID

P04275

Entrez


7450

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1.0 µg/ml, IP: 0.5-1 µg/500 µg protein Lysate, IHC-P: 0.5-1.0 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of human Pa...



Immunohistochemical staining of human To...