

Product Datasheet

CD44 antibody (orb389297)

Description

Mouse monoclonal antibody to CD44

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

Recombinant fragment corresponding to the v9 domain of human CD44

Target

CD44

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

CD44v9 (Marker of Tumor Metastasis)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

80-95kDa

Uniprot ID

P16070

Entrez



960

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human Ce...



Immunohistochemical staining of human To...



Immunohistochemical staining of human Pr...