



Product Datasheet

CD137 antibody (orb389316)


Description	CD137 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD137. Please contact us for the exact immunogen sequence. The peptide is available as orb13629.
Target	CD137
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	28 kDa
Uniprot ID	P20334
NCBI	001552
Dilution Range	WB: 1:200-1500, IHC-P: 1:50-300, IF/ICC: 1:50-300
Expiration Date	12 months from date of receipt.



IHC-P image of rat lymph node tissue using CD137 antibody



Immunofluorescence analysis of rat brain tissue using CD137 antibody



Immunohistochemical staining of paraffin-embedded tissue using CD137 antibody