



Product Datasheet

NOTCH3 antibody (orb389328)


Description	NOTCH3 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human NOTCH3. Please contact us for the exact immunogen sequence. The peptide is available as orb13850.
Target	NOTCH3
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	243 kDa
Uniprot ID	Q9R172 , Q9UM47 , Q61982
NCBI	000426
Dilution Range	WB: 1:100-500, IHC-P: 1:50-500, IF/ICC: 1:50-500
Expiration Date	12 months from date of receipt.



IF analysis of rat lipocyte tissue using...



Immunohistochemical staining of rat lipo...



Western blot analysis of rat lipocyte (l...