



Product Datasheet

GUCY1A2 antibody (orb389331)


Description	GUCY1A2 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GUCY1A2. Please contact us for the exact immunogen sequence. The peptide is available as orb13867.
Target	GUCY1A2
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	82 kDa
Uniprot ID	Q9WV14 , F8VQK3
NCBI	000846
Dilution Range	WB: 1:100-800, IHC-P: 1:100-500, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



IF image of rat ovary tissue using GUCY1...



IF analysis of rat ovary tissue using an...



Immunohistochemical staining of rat pros...